




Accessories for 2.867.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.867.0010

Qty.	Order no.	Description
1 pcs	1.867.0010	<p>867 pH Module</p> <p>Module for pH/ion measurement as supplement to a Titrando or "stand-alone" in combination with a 900 Touch Control.</p> <p>In addition to measurements of pH, temperature, mV, Ipol, Upol and concentration, the pH Module can perform standard additions (manual, dos, autos) and Liquid Handling (add, prep, empty). It enables the use of both conventional and intelligent sensors for measuring. Also integrated in the software is an automatic GLP-compliant pH electrode test.</p> <p>The pH Module has two USB interfaces for connecting printers, barcode readers or sample changers and four MSB interfaces for stirrers or Dosinos (for the addition of auxiliary solutions or for standard addition).</p> <p>Both in conjunction with the Touch Control and also integrated in <i>tiamo</i>[™] full (from 2.0), it is in compliance with GLP and FDA 21 CFR part 11 requirements.</p>
		
1 pcs	6.2103.130	<p>Adapter red 2 mm plug / 4 mm socket</p> <p>For connecting plug B (4 mm) to 2 mm socket.</p>
		
1 pcs	6.2103.140	<p>Adapter black 2 mm plug / B socket 4 mm</p> <p>For connecting plug B (4 mm) to 2 mm socket.</p>
		

Optional accessories

Order no.	Description
-----------	-------------

2.801.0010	801 Stirrer
-------------------	--------------------

Magnetic stirrer without stand for supplementing the Titrino plus, Dosimat plus, Titrandos, Sample Processors, 805 Dosimats and 780/781 pH Meters. With permanently attached cable for MSB (Metrohm Serial Bus).



2.802.0040	802 Stirrer (propeller stirrer) for 804 Ti Stand
-------------------	---

Rod stirrer with 6.1909.010 Propeller stirrer.



2.804.0010	804 Ti Stand without stand rod
-------------------	---------------------------------------

Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.



6.2065.000	Stacking frame for 846 Dosing Interface, 856 Conductivity Module, 867 pH Module
-------------------	--

For fixing the Reagent Organizer on top of the Dosing Interface



6.2151.000 Cable USB A – mini-DIN 8-pin

Controller cable

