








## Accessories for 2.1008.2010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.1008.2010

Qty.	Order no.	Description	
1 pcs	1.1008.0010	<b>Eco Titrator</b>	
		<p>The compact Eco Titrator with integrated magnetic stirrer and touch-sensitive User Interface is ideal for routine analysis. It provides GLP-compliant results with minimum space requirements at all times (approx. DIN A4).</p> <p>Universally compatible with almost all potentiometric titrations, such as, for example, for</p> <ul style="list-style-type: none"> <li>• Food products: Acid content, chloride, Vitamin C, iodine and peroxide number in fats</li> <li>• Water analysis: Carbonate and Ca/Mg hardness, chloride, sulfate, permanganate index</li> <li>• Petrochemistry: Acid/base number, sulfide &amp; mercaptans, chloride, bromine number</li> <li>• Electroplating: Total acid, metal content, chloride</li> <li>• Surfactant analysis: Anionic, cationic and non-ionic surfactants</li> <li>• Photometry with the Optrode: p and m value, metals, water hardness</li> </ul>	
1 pcs	6.00430.100	<b>Ag Titrode</b>	
		<p>Combined silver ring electrode with a pH glass membrane as reference electrode. This maintenance-free electrode is suitable for precipitation titrations when the pH value remains constant (titrant silver nitrate), for example of:</p> <ul style="list-style-type: none"> <li>• chloride, bromide, iodide</li> <li>• sulfides</li> <li>• hydrogen sulfide</li> <li>• mercaptans</li> <li>• cyanides</li> </ul> <p>This electrode is stored in distilled water.</p> <p>Depending on the application, we recommend using an Ag Titrode with a coating (AgBr, AgCl or Ag<sub>2</sub>S coating), which can be ordered accordingly.</p>	
1 pcs	6.03002.210	<b>10 mL cylinder unit for Eco Dosimat / Eco Titrator</b>	
		<p>10 mL cylinder unit for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.</p>	
1 pcs	6.1446.030	<b>Ball stopper / B-14/(15)</b>	
		<p>For the movable holding of buret tips in titration vessels</p>	

<b>1 pcs</b>	<b>6.1446.080</b>	<b>Threaded stopper / M8</b>	
<b>1 pcs</b>	<b>6.1543.200</b>	<b>Titration tip / M6 thread</b> Titration tip with antidiffusion valve.	
<b>1 pcs</b>	<b>6.1546.040</b>	<b>Piston tool for Eco Dosimat/Eco Titrator</b> Tool for removing and mounting the piston.	
<b>1 pcs</b>	<b>6.1602.105</b>	<b>Bottle attachment / GL 45 / green</b> For direct dosing from reagent bottles with GL 45 thread.	
<b>1 pcs</b>	<b>6.1608.023</b>	<b>Amber glass bottle / 1000 mL / GL 45</b> For Exchange Units. Bottle for auxiliary solutions.	

---

**1 pcs**     **6.1619.010**     **Adsorber complete, for Exchange Unit**



---

**1 pcs**     **6.1805.090**     **FEP tubing / M6 / 31 cm**  
With light and kink protection



---

**1 pcs**     **6.1805.100**     **FEP tubing / M6 / 40 cm**  
With light and kink protection



---

**1 pcs**     **6.1819.020**     **Tube**  
For M6 threaded holes. For 6.1602.1XX bottle attachments



---

**1 pcs**    **6.1903.020**    **Stirring bar / 16 mm**  
Stirring Bar with magnetic core, PTFE coated, length 16 mm



---

**1 pcs**    **6.1903.030**    **Stirring bar / 25 mm**  
Stirring bar with magnetic core, PTFE covering, length 25 mm.



---

**1 pcs**    **6.2013.010**    **Clamping ring**  
For support rods with a diameter of 10 mm.



---

**1 pcs**    **6.2016.070**    **Support rod / 400 mm**



---

**1 pcs**    **6.2021.020**    **Electrode holder**  
Electrode holder for 4 electrodes and 2 buret tips



---

**1 pcs**     **6.2023.020**     **Clip for SGJ 14/15**  
Clip for SGJ 14/16



---

**1 pcs**     **6.2104.020**     **Electrode cable / 1 m / F**  
For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



---

**1 pcs**     **6.2164.010**     **Power supply unit 100 - 240 V/24 V DC**  
Power supply unit for 915 KF Ti-Touch and 916 Ti-Touch



---

**1 pcs**     **6.2709.070**     **Guiding sleeve**  
For fixation of buret and dosing tips in titration heads and electrode holders



---

**1 pcs**     **6.6081.010**     **USB flash drive for Eco Titrator**  
Storage medium for Eco Titrator



## Optional accessories

Order no.	Description
<b>2.142.0100</b>	<b>Custom Q3X thermal printer</b>

Compact printer with USB interface for



- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter (cable 6.2151.140)
- Eco Dosimat / Titrator
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer

Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

---