

Accessories for 2.875.9050

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.875.9050

Qty.	Order no.	Description
1 pcs	1.875.9050	875 KF Gas Analyzer without monitor



The 875 KF Gas Analyzer is a completely configured analysis system. The instrument contains a computer with **tiamo™** full titration software installed in addition to an 851 Titrando for automatic coulometric water content determination. Predefined methods for water content determination in liquefied and permanent gases facilitate work with the 875 KF Gas Analyzer.

In addition to the determination of water content in gases, the system can also be used for the determination of water content in liquid and solid samples.

A monitor with keyboard and mouse is required in addition for the control of the system.



Optional accessories

Order no.	Description	
6.7209.000	KF Gas Analyzer equipment for reagent replacement and adding methanol	
	KF Gas Analyzer accessories for automatic reagent replacement or the addition of methanol.	
6.7209.010	KF Gas Analyzer equipment for rinsing with solvent	
	KF Gas Analyzer accessories for rinsing the lines with solvent.	