




Tartozékok: 2.885.0010

Alább a tartozékokat Szállítási csomagra és Opcionális tartozékokra csoportosítva találja. Kérjük, tartsa kéznél ezt a dokumentumot a tartalék anyagok rendeléséhez. Ezek a listák időnként megváltozhatnak.

Szállítási csomag 2.885.0010

Menny.	Rendelész.	Leírás	
1 db	1.885.0010	885 Compact Oven Sample Changer The 885 Compact Oven Sample Changer is used for thermal sample preparation in Karl Fischer titration. In the oven the samples can be heated up to 250°C. The moisture contained in the sample evaporates and is transported by a dry carrier gas into the titration cell, where the analysis takes place.	
1 db	6.1446.170	Heating tube stopper for coulometric KF vessels in connection with heatable tubing connection	
1 db	6.1448.067	Screw cap with septum, 100 pcs. Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces	
2 db	6.1602.145	Drying flask cap Screw cap with GL 45 thread for the 6.1608.050 drying flask.	
2 db	6.1608.050	Drying bottle / 100 mL / GL 45	

1 db **6.1805.010** **FEP tubing / M6 / 13 cm**
With light and kink protection.



1 db **6.1805.060** **FEP tubing / M6 / 60 cm**
With light and kink protection



2 db **6.1805.080** **FEP tubing / M6 / 25 cm**
With light and kink protection.



1 db **6.1808.040** **Thread adapter M6 outer / M8 inner**
Outer thread M6, inner thread M8.



1 db **6.1808.050** **Tubing adapter olive / M8 outer**
1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



2 db **6.1821.040** **Filter tube**
Filter tube for 6.1608.050 Drying bottle. For Rancimats and Karl Fischer ovens.



1 db **6.1830.030** **Heating tubing**
Heating tubing with M6 thread for 860 KF Thermoprep, 885 Compact Oven Sample Processor and 874 USB Oven Sample Processor.



2 db **6.2043.005** **Holding clip for bottles**
Holding spring for reagent bottles in exchange units.



1 db **6.2049.050** **Needle holder with Luer lock**
Needle holder (73 mm) with Luer lock for the Karl Fischer ovens 860, 874 and 885.



1 db 6.2151.100 Adapter USB MINI (OTG) - USB A

For connecting USB instruments.



1 db 6.2420.007 Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces



1 db 6.2621.140 Hexagon key 2.5 mm

1 db 6.2724.010 Dust filter

Dust filter for Rancimats and Karl Fischer ovens.



1 db 6.2739.000 Wrench

For tightening connectors



1 db 6.2811.000 Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



1 db 6.2816.070 Needle for gas inlet

Used with KF Thermoprep.



1 db 6.2816.080 Outlet needle

Used with KF Thermoprep and Oven Sample Processor.

