





## Accessories for 6.1451.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description	
<b>6.1429.040</b>	<b>Reaction vessel short for stability measurements</b>	
	Reaction vessel for stability measurements with Rancimat and PVC Thermomat. Set of 117 pieces.	
<b>6.1429.050</b>	<b>Reaction vessel long for stability measurements</b>	
	Reaction vessel for stability measurements involving volatile oils, mineral oils and biodiesel blends with the Biodiesel Rancimat. In combination with 6.2418.130. Set of 100 pieces.	
<b>6.2418.100</b>	<b>Air tube for oil measurements</b>	
	Set of 117 pieces. In combination with 6.1429.040.	
<b>6.2418.130</b>	<b>Air tube long for biodiesel measurements</b>	
	Set of 100 pieces. In combination with 6.1429.050.	