





Accessories for 2.909.0014

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.909.0014

Qty.	Order no.	Description	
1 pcs	1.909.0014	909 UV Digester (230 V)	
		Digestion instrument for UV photolysis of water samples with low to medium organic load. For sample preparation in trace element determination by means of voltammetry, ion chromatography and spectroscopy (AAS, ICP). Integrated instrument with operating unit and wet end. With air cooling and automatic control of digestion temperature and time. For 12 samples with a maximum of 12 mL sample volume each. Instrument for 220 - 240 V and 50 - 60 Hz.	
1 pcs	6.1110.010	Pt1000 temperature sensor for 909 UV Digester	
		Pt1000 temperature sensor (class B) made of glass with fixed cable (0.5 m) for connecting to a 909 UV Digester.	
12 pcs	6.1446.100	PTFE stopper for 12 mL UV quartz sample vessel	
		For 6.2414.000 quartz sample vessel	
1 pcs	6.1446.190	Thermometer stopper for Pt1000 temperature sensor	
		For 6.1110.010 Pt1000 temperature sensor and 6.2414.000 12 mL quartz sample vessel.	

1 pcs **6.2008.040** **Storage vessel**

Together with 6.2008.050 storage vessel holder provides a support for the electrode on 6.3032.XXX dosing units.



1 pcs **6.2041.240** **Sample holder for 12 quartz sample vessels 12 mL**

For 909 UV Digester.



1 pcs **6.2225.000** **UV warning plate**

For 909 UV Digester



12 pcs **6.2414.000** **Quartz sample vessel 12 mL**

For UV Digester



1 pcs **6.2621.100** **3 mm hex key for IC Sample Processors**



1 pcs **6.2745.100** **UV protective shield**

For 909 UV Digester.



1 pcs **6.2804.090** **UV mercury vapor lamp**

Medium pressure mercury vapor lamp for 909 UV Digester.



Optional accessories

Order no.	Description
6.2042.050	Spacer for Pt1000 temperature sensor 