

## Accessories for BWT-840000972

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery BWT-840000972

Qty.	Order no.	Description	
1 pcs	<b>BWT-840000101</b>	<b>Laser safety goggles (785 nm)</b>	
		Laser safety goggles for 785 nm lasers – suitable only for Class IIIB laser systems from B&W Tek.	
1 pcs	<b>BWT-840000295</b>	<b>Vials for liquids, 15 mm in diameter</b>	
		Borosilicate glass vials (diameter 15 mm), package of 6, with pouch.	
1 pcs	<b>BWT-840000329</b>	<b>Pouch for 15 mm Vials</b>	
		Pouch for 15 mm vials.	
1 pcs	<b>BWT-840000403</b>	<b>10.8 V rechargeable lithium ion battery</b>	
		Rechargeable lithium-ion battery, for use with NanoRam-1064, TacticID-GP Plus, TacticID-N Plus, TacticID-1064.	
1 pcs	<b>BWT-840000406</b>	<b>TacticID Vial-Holder Adapter</b>	
		TacticID Vial-Holder Adapter, compatible with 15 mm-diameter vials.	
1 pcs	<b>BWT-840000414</b>	<b>18 V Power Adapter</b>	
		Power adapter for 18 V DC, 1.67 A. Plugs for USA, EU, UK and AU included in the scope of delivery.	

---

**1 pcs**    **BWT-840000416**    **TacticID Point-and-Shoot adapter**

TacticID Point-and-Shoot adapter.



---

**1 pcs**    **BWT-840000417**    **TacticID Validation Cap**

TacticID Validation Cap with integrated ASTM-standard polystyrene. For use with TacticID-N Plus and TacticID-GP Plus.



---

**1 pcs**    **BWT-840000418**    **TacticID-GP Protection Cap**

TacticID-GP shaft protective cap.



---

**1 pcs**    **BWT-840000422**    **TacticID Carrying Case**

TacticID hard plastic rugged carrying case.



---

**1 pcs**    **BWT-840000491**    **TacticID Right-Angle Adapter**

TacticID Right-Angle Adapter.



---

**1 pcs**    **BWT-840000858**    **USB cable**

USB cable supporting micro USB port up to type A USB 2.0.







---

**1 pcs**    **BWT-840000861**    **Micro USB flash drive**

TacticID micro USB flash drive.



1 pcs	<b>BWT-840000912</b>	<b>O-ring replacement kit for TacticID instruments</b>	
<p>TacticID O-ring replacement kit – includes: O-rings, ID 9.00 mm x15; O-Rings, ID 10.00 mm x15; O-rings, ID 13.00 mm x10 for the shaft in handheld Raman instruments.</p>			
1 pcs	<b>BWT-840000948</b>	<b>USB OTG Host Cable</b>	
<p>USB OTG host cable for data transmission between a USB flash drive and a handheld system.</p>			
1 pcs	<b>BWT-840000974</b>	<b>Complete Raman Spectral Library (TacticID Plus)</b>	
<p>Complete Raman Spectral Library for TacticID Plus with &gt; 10,000 spectra of common chemicals, narcotics, solvents and plastics. Includes NFPA704 and GHS safety instructions.</p>			
1 pcs	<b>BWT-840000975</b>	<b>Common Chemicals, Narcotics, Explosives, Solvents and Plastics Spectral Library (TacticID Plus)</b>	
<p>Raman Spectral Library for TacticID Plus with more than 1,200 spectra of common chemicals, narcotics, explosives, solvents and plastics. Includes NFPA704 and GHS safety instructions. May be subject to export control.</p>			
1 pcs	<b>BWT-840000980</b>	<b>TacticID Plus Manuals</b>	
<p>TacticID Plus Product Manual Set including Quick Reference Guide and TID Installation Instructions.</p>			
1 pcs	<b>BWT-840000981</b>	<b>TacticID Plus Software</b>	
<p>USB flash drive containing TacticID Plus TID software for PC, TID installation package and TacticID Plus Product Manuals.</p>			

## Optional accessories

Order no.	Description	
<b>BWT-840000294</b>	<b>Battery Charging Cradle</b>  Lithium-ion battery charging cradle.(for charging NR-BAT and TID-BAT only).	
<b>BWT-840000353</b>	<b>Shaft Protection Cap</b>  Handheld Raman shaft protective cap, package of 10.	
<b>BWT-840000361</b>	<b>Handheld Laser Power Meter</b>  Laser Power Meter, Pen Type, Calibrated. Range: 5 micro-Watts to 1 Watt. Used in OQ/PQ and annual certification laser power measurement.	
<b>BWT-840000428</b>	<b>Wireless Router</b>  Wireless, battery-powered portable router; 3G/4G connectivity.	
<b>BWT-840000498</b>	<b>Disposable Battery Cartridge</b>  TactiID battery cartridge for use with six type 123A disposable batteries (e.g. Surefire). Batteries not included.	
<b>BWT-840000524</b>	<b>TacPac™ Adapter</b>  SERS analysis adapter for use with TacPac™-P SERS substrates.	
<b>BWT-840000525</b>	<b>TacPac™ SERS substrate (package of 10)</b>  Performance-grade, flexible and sensitive SERS substrates in strip form for dip or liquid pipette sample techniques. Package with 10 units.	
<b>BWT-840000526</b>	<b>TacPac™ SERS substrate (package of 30)</b>  Performance-grade, flexible and sensitive SERS substrates in strip form for dip or liquid pipette sample techniques. Package with 30 units.	

---

**BWT-840000527 TacPac™ sample preparation kit (package of 10)**

Disposable plastic single-use laboratory spatula, pipettes and PP graduated microcentrifuge tubes (2 ml) used to prepare samples for testing with TacPac™ SERS. Package of 10.



---

**BWT-840000528 TacPac™ Kit for TacticID**

SERS analysis kit for use with Raman handheld spectrometers from B&W Tek. For Raman signal amplification and fluorescence reduction of hard-to-measure samples. For use with 785 nm laser excitation. Includes: tweezers, TacPac adapter, 10 SERS strips.



---

**BWT-840000585 TacPac™ sample preparation kit (package of 30)**

Disposable plastic single-use laboratory spatula, pipettes and PP graduated microcentrifuge tubes (2 ml) used to prepare samples for testing with TacPac™ SERS. Package of 30.



---

**BWT-840000603 Pouch for TacticID Adapters**

TacticID (785 nm) Pouch for Adapters



---

**BWT-840000620 TacPac™ SERS substrate (package of 100)**

Performance-grade, flexible and sensitive SERS substrates in strip form for dip or liquid pipette sample techniques. Package with 100 units.



---

**BWT-840000621 TacPac™ sample preparation kit (package of 100)**

Disposable plastic single-use laboratory spatula, pipettes and PP graduated microcentrifuge tubes (2 ml) used to prepare samples for testing with TacPac™ SERS. Package of 100.



---

**BWT-840000669 TacticID Car Charger**

TacticID mobile car charger.



---

**BWT-840000938 MicroUSB-USB Flash Drive**

Micro USB to USB flash drive for data storage.



---

**BWT-840000976 Common Chemicals, Narcotics, Solvents and Plastics Spectral Library (TacticID Plus)**

Raman Spectral Library for TacticID Plus with more than 1,180 spectra of common chemicals, narcotics, solvents and plastics. Includes NFPA704 and GHS safety instructions.



---

**BWT-840000977 Raman Explosives Spectral Library (TacticID Plus)**

Raman Explosives Spectral Library TacticID Plus including > 180 spectra. Includes NFPA704 and GHS safety instructions.

Note: This item may be subject to Export control to some territories.



---

**BWT-840000978 Raman Gemstone Spectral Library (TacticID Plus)**

Raman Gemstone Spectral Library for TacticID Plus including 467 spectra.



---

**BWT-840000979 Raman Spectral Library of Chemical Warfare Agents (TacticID Plus)**

Raman Spectral Library of Chemical Warfare Agents for TacticID Plus. Includes NFPA704 and GHS safety instructions. Note: This item may be subject to Export control to some territories.



---

**BWT-840000982 TacticID Plus Remote Operation**

Tablet with Android operating system pre-loaded with the TOS Plus Remote application installed for remote data acquisition and access. Includes camera and GPS functions. TID-WPR, Wireless router included.

