






Accessories for A629280110

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery A629280110

Qty.	Order no.	Description	
1 pcs	6.7440.480	NIRS Process, I/O Communication box I/O Communication box for process instruments with 24 analog and 24 digital channels	
1 pcs	6.7450.070	NIRS 99% reflection and wavelength standard process, 1 inch diameter NIRS 99% reflection and wavelength standard for process instrument with 1 inch probe diameter	
1 pcs	6.7450.080	NIRS 99% reflection and wavelength standard for process instrument, 0.5 inch probe diameter NIRS 99% reflection and wavelength standard for process instrument, 0.5 inch probe diameter	
1 pcs	6.7450.090	NIRS 99% reflection and wavelength standard for process instrument, 1 inch diameter, 45° bevel NIRS 99% reflection and wavelength standard for process instrument with 1 inch diameter and 45° angle	

Optional accessories

Order no.	Description	
A649281190	NIRS XDS Process Contact Lamp 50W, 12 Volt- White lead	
Replacement lamp with contact for NIRS Process instruments. 50W, 12 Volt- White lead.		
