

Accessories for BWT-84000920

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery BWT-84000920

Qty.	Order no.	Description	
1 pcs	BWT-84000102	Laser safety goggles (1064 nm)	
		Laser safety goggles for lasers with excitation at 1,064 nm – suitable for Class IIIB laser systems from B&W Tek.	
1 pcs	BWT-84000295	Vials for liquids, 15 mm in diameter	
		Borosilicate glass vials (diameter 15 mm), package of 6, with pouch.	
1 pcs	BWT-84000329	Pouch for 15 mm Vials	
		Pouch for 15 mm vials.	
1 pcs	BWT-84000342	Right-Angle Adapter	
		Right-angle adapter for hands-free measurement using handheld Raman instruments. Compatible with NanoRam, NanoRam-1064 and TacticID-1064 models.	
1 pcs	BWT-84000343	Point-and-Shoot-Adapter	
		Point-and-shoot adapter for handheld Raman instruments. Compatible with NanoRam, NanoRam-1064 and TacticID-1064 models.	
1 pcs	BWT-84000344	Vial Holder Adapter	
		Vial holder adapter, compatible with 15 mm-diameter vials.	
1 pcs	BWT-84000345	Handheld Raman bottle adapter	
		Handheld Raman accessories: bottle adapter	

1 pcs	BWT-84000348	Polystyrene Validation Cap	<p>Polystyrene Validation cap with integrated ASTM standard polystyrene.</p> 
1 pcs	BWT-84000354	Pouch for adapters	<p>Pouch for NanoRam adapters.</p> 
1 pcs	BWT-84000403	10.8 V rechargeable lithium ion battery	<p>Rechargeable lithium-ion battery, for use with NanoRam-1064, TacticID-GP Plus, TacticID-N Plus, TacticID-1064.</p> 
1 pcs	BWT-84000414	18 V Power Adapter	<p>Power adapter for 18 V DC, 1.67 A. Plugs for USA, EU, UK and AU included in the scope of delivery.</p> 
1 pcs	BWT-84000419	TacticID Stylus/pen combo	<p>Stylus/pen combo for touch screens (package of 5).</p> 
1 pcs	BWT-84000858	USB cable	<p>USB cable supporting micro USB port up to type A USB 2.0.</p> 
1 pcs	BWT-84000911	TacticID Large Carrying Case	<p>Carrying case for TacticID-1064 and TacticID ST.</p> 

1 pcs BWT-84000912 O-ring replacement kit for TacticID instruments

TacticID O-ring replacement kit – includes: O-rings, ID 9.00 mm x15; O-Rings, ID 10.00 mm x15; O-rings, ID 13.00 mm x10 for the shaft in handheld Raman instruments.



1 pcs BWT-84000939 TID-EX Software Drive and Product Manual

USB flash drive containing TacticID-1064 PC software TID-EX installation package and TacticID-1064 product manuals.



1 pcs BWT-84000948 USB OTG Host Cable

USB OTG host cable for data transmission between a USB flash drive and a handheld system.



Optional accessories

Order no.	Description	
BWT-840000294	Battery Charging Cradle	
	Lithium-ion battery charging cradle.(for charging NR-BAT and TID-BAT only).	
BWT-840000353	Shaft Protection Cap	
	Handheld Raman shaft protective cap, package of 10.	
BWT-840000361	Handheld Laser Power Meter	
	Laser Power Meter, Pen Type, Calibrated. Range: 5 micro-Watts to 1 Watt. Used in OQ/PQ and annual certification laser power measurement.	
BWT-840000428	Wireless Router	
	Wireless, battery-powered portable router; 3G/4G connectivity.	
BWT-840000498	Disposable Battery Cartridge	
	TacticID battery cartridge for use with six type 123A disposable batteries (e.g. Surefire). Batteries not included.	
BWT-840000669	TacticID Car Charger	
	TacticID mobile car charger.	
BWT-840000969	Contact Immersion Probe (1064 nm)	
	Handheld Raman Spectrometer 1064 accessories: Contact immersion probe, 6" length, working distance approx. 1 mm, operating pressure up to 2 bar at room temperature.	

BWT-84000970

Handheld Raman large bottle adapter (1064)

Handheld Raman Spectrometer 1064 accessories: Large bottle adapter with 10 mm working distance.

